

Selectra - Flexor Systems

# ProActive

MAINTENANCE KITS



from  
ELITech  
Clinical Systems



Vital Scientific BV - PO Box 100 - 6950 AC - DIEREN - The Netherlands

Tel: +31 (0)313-430 500 - Fax: +31 (0)313-427 807 - E-mail: [vital.service@elitechgroup.com](mailto:vital.service@elitechgroup.com) - Website: [www.elitechgroup.com](http://www.elitechgroup.com)

*ELITech helps your laboratory take the guesswork out of periodic maintenance to keep your instrument running at peak performance.*

The ProActive Maintenance Kits contain all the necessary replacement parts for a specific scheduled maintenance. The content of each kit is carefully defined to match to Selectra and Flexor Systems bi annual and annual maintenance schedule.

Each kit is specified based on ELITech's extensive system and maintenance experience and is designed to be cost effective.

- ✓ Pre-defined 1 year single kit. Genuine factory service parts provided according to the maintenance schedule
- ✓ Kit prices are more economical than the purchase prices of individual parts and designed to reduced maintenance costs
- ✓ Easy-to-plan annual maintenance parts budget
- ✓ Minimal spare parts ordering and inventory management
- ✓ Effective preventive maintenance improves system reliability and enhances system uptime
- ✓ Up to date preventative maintenance instructions guide customers to replace the right parts
- ✓ Every kit has catalogue number that makes ordering straightforward and easy.

*Contact your local distributor or sales representative for further details*

Kind regards,  
The System Support Group

## Contents:

### Introduction

<b>Microlab Series</b>	1
6002-469: Yearly Maintenance Kit	
<b>Selectra ProXS – Flexor EL80</b>	2
6003-253: Yearly Maintenance Kit	
<b>Junior Series</b>	3
6003-475: Yearly maintenance Kit	
6002-898: Supplementary Two Yearly Maintenance Kit	
<b>Selectra ProS – Flexor EI150</b>	4
6003-539: Yearly Maintenance Kit	
6002-898: Supplementary Two Yearly Maintenance Kit	
<b>E Series</b>	5
6003-476: Yearly Maintenance Kit	
6002-152: Supplementary Two Yearly Maintenance Kit	
<b>Selectra ProM – Flexor EL200</b>	6
6003-383: Yearly Maintenance Kit	
6002-152: Supplementary Two Yearly Maintenance Kit	
<b>XL Series</b>	7
6003-447: Yearly Maintenance Kit	
6002-680: Supplementary Two Yearly Maintenance Kit	



**6002-469: Yearly Maintenance Kit**

Part number	Description	Qty
1513-009	Tubing pump unit (need 2x75mm),	1
3029-039	Pump O-ring	1
3380-018	Lamp, clear quartz 12V/20W	1
4100-347	Bellows	1
4104-134	Sipper tubing (255mm)	1
6002-330	Flowcell tube assy	1
6002-331	Waste tube assy	1



SELECTRA PRO XS  
FLEXOR EL 80

### 6003-253: Yearly Maintenance Kit

Part number	Description	Qty
1513-001	Tube, PVC 2x4 mm clear 250mm	1
3053-020	Ball bearing 8x16x5mm	1
3064-041	Drivingbelt mixer	2*
3066-083	Water filter	2*
3066-155	Syringe Long Life 100 uL	1
3066-156	Syringe Long Life 1 mL	1
3380-018	Lamp, clear quartz 12V/20W	1
3913-058	Service kit water/cooling pump	3**
6002-056	Service kit water pump	2**
6002-776	Needle tube assy	1

\* Amount is sufficient for multiple user/service maintenance for a period of one year

\*\* The number to be used is depending on the configuration of the instrument because it is possible to have EPDM or NON-EPDM valves in the water pump. By the end of 2012 we will stop adding both type of membranes. Only the EPDM membranes will be provided by the end of 2012.



**6003-475 Yearly Maintenance Kit Junior V2**

Part number	Description	Qty.
1513-001	Tube, PVC 2x4 mm clear 250mm	1
3053-020	Ball bearing 8x16x5mm	1
3064-041	Drivingbelt mixer	2*
3066-083	Water filter	2*
3066-155	Syringe Long Life 100 uL	1
3066-156	Syringe Long Life 1 mL	1
3380-018	Lamp, clear quartz 12V/20W	1
3913-058	Service kit water/cooling pump	2**
6001-963	Filter block assy	4*
6002-056	Service kit water pump	2**
6002-775	Connection tube assy	1
6002-776	Needle tube assy	1
6002-777	Rinse tube assy	1
6003-153	Membrane kit EPDM (Vacuum Pump)	1

\* Amount is sufficient for multiple user/service maintenance for a period of one year

\*\* The number to be used is depending on the configuration of the instrument because it is possible to have EPDM or NON-EPDM valves in the water pump. By the end of 2012 we will stop adding both type of membranes. Only the EPDM membranes will be provided by the end of 2012.

**6002-898: Supplementary Two Yearly Maintenance kit**

Part number	Description	Qty.
3066-072	Valve block	1
6003-154	Chamber kit EPDM (Vacuum Pump)	2




### 6003-539 Yearly Maintenance Kit

Part number	Description	Qty.
1513-001	Tube, PVC 2x4 clear 250mm	1
3053-020	Ball bearing 8x16x5mm (mixer)	1
3064-041	Drivingbelt mixer	2*
3066-083	Water filter	2*
3066-155	Syringe Long Life 100 uL	1
3066-156	Syringe Long Life 1 mL	1
3380-018	Lamp, clear quartz 12V/20W	1
3913-058	Service kit water/cooling pump	2*
6001-963	Filter block assy	4*
6002-775	Connection tube assy	1
6002-776	Needle tube assy	1
6002-777	Rinse tube assy	1
6003-153	Membrane kit EPDM (Vacuum Pump)	1

\* Amount is sufficient for multiple user/service maintenance for a period of one year

### 6002-898: Supplementary Two Yearly Maintenance Kit

Part number	Description	Qty.
3066-072	Valve block	1
6003-154	Chamber kit EPDM (Vacuum Pump)	2


**6003-476: Yearly Maintenance Kit E series V2**

Part number	Description	Qty.
3053-020	Ball bearing 8x16x5mm (sample mixer)	1
3053-038	Ball bearing 7x14x5mm (reagent mixer)	1
3064-041	Drivingbelt mixer	4*
3066-083	Water filter	2*
3066-155	Syringe Long Life 100 uL	1
3066-156	Syringe Long Life 1 mL	1
3370-730	Code marker "R"	2
3370-731	Code marker "S"	2
3380-018	Lamp, clear quartz 12V/20W	1
3913-058	Service kit water/cooling pump	1**
6001-405	Reagent tube assy	1
6001-428	Sample tube assy	1
6001-431	Tube assy dispenser	2
6001-861	Tube assy pump unit	1
6001-963	Filter block assy	4*
6002-056	Service kit water pump	1**
6003-153	Membrane kit EPDM (Vacuum Pump)	1

\* Amount is sufficient for multiple user/service maintenance for a period of one year

\*\* The number to be used is depending on the configuration of the instrument because it is possible to have EPDM or NON-EPDM valves in the water pump. By the end of 2012 we will stop adding both type of membranes. Only the EPDM membranes will be provided by the end of 2012.

**6002-152: Supplementary Two Yearly Maintenance Kit**

Part number	Description	Qty.
3066-072	Valve block	2
6003-154	Chamber kit EPDM (Vacuum Pump)	2





**6003-383: Yearly Maintenance Kit**

Part number	Description	Qty
3053-020	Ball bearing 8x16x5mm (sample mixer)	1
3053-038	Ball bearing 7x14x5mm (reagent mixer)	1
3064-041	Drivingbelt mixer	4*
3066-083	Water filter	2*
3066-155	Syringe Long Life 100 uL	1
3066-156	Syringe Long Life 1 mL	1
3380-018	Lamp, clear quartz 12V/20W	1
3913-058	Service kit water/cooling pump	2*
6001-963	Filter block assy	4*
6003-153	Membrane kit EPDM (Vacuum Pump)	1
6003-336	Reagent tube assy	1
6003-337	Sample tube assy	1

\* Amount is sufficient for multiple user/service maintenance for a period of one year

**6002-152: Supplementary Two Yearly Maintenance Kit**

Part number	Description	Qty
3066-072	Valve block	2
6003-154	Chamber kit EPDM (Vacuum Pump)	2



**XL series**

**6003-477: Yearly Maintenance Kit**

Part number	Description	Qty.
3053-020	Ball bearing 8x16x5mm (sample mixer)	2
3053-038	Ball bearing 7x14x5mm (reagent mixer)	2
3064-041	Drivingbelt mixer	8*
3066-083	Water filter	2*
3066-155	Syringe Long Life 100 uL	2
3066-156	Syringe Long Life 1 mL	2
3370-764	Dust filter	2
3380-018	Lamp, clear quartz 12V/20W	2
3913-058	Service kit water/cooling pump	1**
6001-963	Filter block assy	8*
6002-056	Service kit water pump	1**
6002-663	Sample tube assy	2
6002-664	Reagent tube assy	2
6003-153	Membrane kit EPDM (Vacuum Pump)	2

\* Amount is sufficient for multiple user/service maintenance for a period of one year

\*\* The number to be used is depending on the configuration of the instrument because it is possible to have EPDM or NON-EPDM valves in the water pump. By the end of 2012 we will stop adding both type of membranes. Only the EPDM membranes will be provided by the end of 2012.

**6002-680: Supplementary Two Yearly Maintenance Kit**

Part number	Description	Qty.
3066-072	Valve block	4
6003-154	Chamber kit EPDM (Vacuum Pump)	4